



Innate Immunity to Life!



3 for 2 *

*Promo code 3IAX2 required when ordering!
Valid until April 30, 2014.

Order 3 ODNs listed below &
receive the cheapest for free!

The ODN Experts – Competitive Prices

- Endotoxin-free and sterile
- Potent TLR9 ligands
- Tested on TLR9 KO murine macrophages
- BULK available for *in vivo* studies
- ddWater Endotoxin-free (sterile) included
- Unique TLR9 Explorer Sets™
- Inhibitory ODNs (iODNs) available
- Negative ODN Controls available

PID	PRODUCT NAME		SIZE
ODNs			
IAX-200-005	ODN 2216 (Type A) Endotoxin-free (sterile)	STANDARD	100 µg 1 mg 3 x 1 mg
IAX-200-003	ODN 1585 (Type A) Endotoxin-free (sterile)		100 µg 1 mg 3 x 1 mg
IAX-200-001	ODN 1668 (Type B) Endotoxin-free (sterile)	STANDARD	100 µg 1 mg 3 x 1 mg
IAX-200-002	ODN 1826 (Type B) Endotoxin-free (sterile)	STANDARD	100 µg 1 mg 3 x 1 mg
IAX-200-006	ODN 2006 (Type B) Endotoxin-free (sterile)		100 µg 1 mg 3 x 1 mg
IAX-200-004	ODN M362 (Type C) Endotoxin-free (sterile)		100 µg 1 mg
IAX-200-007	ODN 2395 (Type C) Endotoxin-free (sterile)		100 µg 1 mg 3 x 1 mg
IAX-200-024	AT-ODN-1 Endotoxin-free (sterile)		100 µg
IAX-200-025	AT-ODN-2 Endotoxin-free (sterile)		100 µg
IAX-200-026	AT-ODN-3 Endotoxin-free (sterile)		100 µg
iODNs			
IAX-200-050	iODN (inhibitory ODN) 2088 Endotoxin-free (sterile)	STANDARD	100 µg 1 mg 3 x 1 mg
IAX-200-051	iODN (inhibitory ODN) (ttaggg)4 Endotoxin-free (sterile)		100 µg 1 mg 3 x 1 mg
IAX-200-052	G-type iODN (inhibitory ODN) Endotoxin-free (sterile)		100 µg 1 mg
IAX-200-053	Mini-iODN (inhibitory ODN) Endotoxin-free (sterile)		100 µg 1 mg
IAX-200-054	Mega-iODN (inhibitory ODN) Endotoxin-free (sterile)		100 µg 1 mg
IAX-200-055	Duo-iODN (inhibitory ODN) Endotoxin-free (sterile)		100 µg 1 mg
Controls			
IAX-200-200	ODN 1720 (Control for ODN 1668) Endotoxin-free (sterile)		100 µg
IAX-200-201	ODN 1982 (Control for ODN 1826) Endotoxin-free (sterile)		100 µg
IAX-200-202	Neutral-ODN (Control for iODNs) Endotoxin-free (sterile)		100 µg 1 mg 3 x 1 mg
IAX-200-203	ODN 2118 (Control for ODN 1585) Endotoxin-free (sterile)		100 µg
IAX-200-204	ODN M383 (Control for ODN M362) Endotoxin-free (sterile)		100 µg
IAX-200-205	ODN 2243 (Control for ODN 2216) Endotoxin-free (sterile)		100 µg
IAX-200-206	ODN 2137 (Control for ODN 2006) Endotoxin-free (sterile)		100 µg